

HEALTH CARE STERILIZERS

- **10 FDA Validated Cycles**
- **4 User Adjustable Cycles**
- **Gravity, Fluid and Vacuum Cycles**
- **Built-in Menus and Intuitive Screens**
- **Sealed Keypad and 25 line CRT Display**
- **Permanent Record – Impact Printer**
- **High Vacuum Exhaust**
- **Built-In Steam Generator or user steam**
- **Automatic Condensate Cooling with Water Savings**
- **ASME Stamped Pressure Vessel(s)**
- **UL, cUL listed. Meets ANSI/AAMI Stds**
- **Industry Recognized Reliability**
- **Built-in Service Routines**
- **Effluent Decontamination Option**



Consolidated Healthcare Sterilizers provide a cost effective and reliable solution for Health Care facilities operating rooms, medical clinics and central sterilization supplies. FDA validated cycles plus user define cycles provide a maximum of flexibility and minimum of operator involvement. All Healthcare sterilizers include a high vacuum system for fast processing and automatic condensate cooling and water savings, plus options for steam sources, transport systems, special enclosures and more.

All Healthcare Sterilizers incorporate a long lasting, solid silicone door gasket that can be changed without tools or adhesives, and non-proprietary valves that are often available from local suppliers. Together with our philosophy of no required routine maintenance, they add up a very low cost of ownership and maximum uptime.

With over 50 years of heritage, these sterilizers feature Consolidated's industry renowned quality which together with our network of local dealers and service companies insure the lowest cost of ownership in the industry while maintaining the highest standards of sterilization.



STANDARD FEATURES

- High Vacuum
- Automatic Condensate Exhaust Cooling and water savings
- Stainless Steel Chamber
- Stainless Steel Cabinet
- Full Length Bottom Shelf and Two Piece Top Shelf
- Emergency Manual Override
- ASME Certified Vessel and Steam Generator,
- UL, cUL listed

OPTIONAL FEATURES

- Built-in Electric Steam Generator
- Chamber Flood Switch
- Effluent Decontamination System
- Automatic Jacket Blowdown
- Automatic Generator Blowdown
- Automatic Feed Water Pump
- Recess Paneling
- Rack with 2 removable shelves
- Loading Cart and Transfer Carriage
- Electro-polish
- Passivation
- Perforated Instrument Trays

CYCLES AND PARAMETERS

Key	Cycle Type	Pressure (PSI)	Temp (F)	Time (min)	Dry Time (min)
1	Gravity-30	20.0	250.0	30.0	30.0
2	Gravity-10	32.0	270.0	10.0	30.0
3	Flash-3	32.0	270.0	3.0	1.0
4	Flash-10	32.0	270.0	10.0	1.0
5	Fluids-30	20.0	250.0	30.0	Exhaust Rate: 1PSI/min
6	Fluids-45	20.0	250.0	45.0	Exhaust Rate: 1PSI/min
7	Prevac-Unwrap	32.0	270.0	4.0	1.0
8	Prevac-Wrap	32.0	270.0	4.0	1.0
9	Bowie-Dick	32.0	270.0	3.0	
10	Vacuum-Leak Test	0.0	0	15.0	
11	Gravity-User	20.0	250.0	20.0	25.0
12	Flash-User	32.0	270.0	3.0	10.0
13	Fluids-User	20.0	250.0	35.0	Exhaust Rate: 1PSI/min
14	Prevac-User	32.0	270.0	4.0	10.0

Cycles 1-10 are fixed per FDA Validation. Cycles 11-14 are user adjustable.