



SterilCARE™

The SterilCARE™ Proactive Maintenance Program (PMP) is designed to help keep your autoclave running at top performance for many years. Choose from four SterilCARE™ plans, or contact us to arrange for a custom configured maintenance program.

Choose the SterilCARE™ plan that best suits your facility

Typical Daily Usage	SterilCARE™ Plan Name (recommended)	SterilCARE™ Plan Components				Consolidated Order Number
		Scheduled Inspection	On-Site Service Call	Parts	Travel and Labor Rate	
1 to 2 cycles	Basic	1 per year			Preferred, Discounted	SCB
3 to 4 cycles	Deluxe	2 per year	1 per year			SCD
5 to 6 cycles	Premium	3 per year	2 per year			SCP
> 6 cycles	Elite	4 per year	Unlimited	Included		SCE
	Custom	Contact Consolidated to develop a custom configured SterilCARE™ plan for your facility.				

SterilCARE™ PartsPAX™ Maintenance Kits

As part of the SterilCARE™ Proactive Maintenance Program (PMP), Consolidated offers the convenience of PartsPAX™ Maintenance Kits for your autoclave. PartsPAX™ packages **reduce costs** and include some of the most commonly needed items to keep your autoclave running smoothly and reliably. Consolidated offers two different levels to fit your needs, **Basic** and **Deluxe**. An Options PartsPAX™ Kit includes an assortment of additional replacement items required on occasion.

For Steam Sterilizers	PMP PartsPAX™	
Consolidated Order Number	14-001	14-002
	Basic	Deluxe
Steam Valve Repair Kit (Chamber)	Yes	Yes
Door Gasket	Yes	Yes
Internal Door Lock Kit (Diaphragm, Spring, Gasket)	Yes	Yes
Safety Valve (Jacket)	Yes	Yes
Safety Valve (Chamber)*	Yes	Yes
Steam Valve Repair Kit (Jacket)*	Yes	Yes
Steam Trap Kit (Chamber)		Yes
Steam Trap Kit (Jacket)*		Yes
Fast Exhaust Valve Repair Kit		Yes
Slow Exhaust Valve Repair Kit		Yes
Jacket Exhaust Valve Repair Kit*		Yes
Water Valve Repair Kit*		Yes

For Steam Generators	PMP PartsPAX™	
Consolidated Order Number	14-004	14-005
	Basic	Deluxe
Generator Gasket	Yes	Yes
Water Fill Valve Repair Kit	Yes	Yes
Contactactor	Yes	Yes
Water Level Control Gasket*	Yes	Yes
Safety Valve	Yes	Yes
Sight Glass with 2 Gaskets	Yes	Yes
Sight Glass Valves (Set of 2)	Yes	Yes
Generator Blowdown Valve Repair Kit*	Yes	Yes
Water Level Controller		Yes
Wires and Buss Bars		Yes
Heating Element Set with Mounting Flange and Gasket		Yes
Over-temperature Controller (Chromalox)*		Yes
Over-temperature Probe (Chromalox)*		Yes

Options PartsPAX™

The Options PartsPAX™ Kit includes an assortment of additional replacement items typically required over the life of the sterilizer.

For Steam Sterilizers	Options PartsPAX™
Consolidated Order Number	14-003
Air Filter	Included
Check Valves (3)	Included
Chamber Drain Strainer	Included
High Temperature Door Grease (2 tubes)	Included
Vacuum Release Valve Repair Kit*	Included

Ordering Information

To order call your area Consolidated Parts Distributor or contact Consolidated. Please provide your sterilizer serial number when ordering. This program does not permit substitutions or deletions.



* Supplied when applicable to unit

